



SEED GERMINATION

Seed germination chambers and rooms are specifically designed to provide optimum germination conditions.

SEED GERMINATION



Model #	Growth Area	Growth Height	Lighting PPFD	Temperature		Exterior Dimensions Width x Depth x Height
				Lights On	Lights Off	
CHAMBERS						
SGC-3	0.3m ² x 5 (3.4ft ² x 5)	230mm x 5 (9" x 5)	18μmol m ⁻² s ⁻¹	10°C to 40°C	4°C to 40°C	735mmW x 915mmD x 2085mmH (29"W x 36"D x 82"H)
ROOMS						
SGR-XXX	Customized to your research needs contact BioChambers					

Note a: Light intensity measurements at the shelf at 25°C.



Biochambers Seed Germination Products version 2023-07A.
Our policy of continuous product improvement will occasionally result in changes to product specifications without notice.
©BIOCHAMBERS INCORPORATED 2023. ALL RIGHTS RESERVED PRINTED IN CANADA

www.biochambers.com